

Real-time Quantification of Impurities in Food-Grade CO₂ with PTR-MS

Jens Herbig, Lukas Märk, Felix Piel, Alfons Jordan, Gernot Hanel, Eugen Hartungen, Stefan Jaksch, Simone Jürschik and Philipp Sulzer

IONICON Analytik GmbH., Eduard-Bodem-Gasse 3, 6020 Innsbruck, AUSTRIA



IONICON

Introduction

Proton-Transfer-Reaction Mass Spectrometry (PTR-MS) is a well-established technology for **real-time analysis** of volatile organic compounds (VOCs), which is characterized by high sensitivity, low Limit-of-Detection (LoD), short response time and direct injection without any sample preparation.

Typical applications can be found in various fields such as **environmental sciences, food and flavor** research and many more [1].

In recent years analytical challenges in trace gas detection have become more and more demanding and the need for instrumental performance has continuously grown.

A powerful option to meet these requirements is the **optimization of the ion transfer system** between the PTR ionization region and the TOF mass spectrometer, which can **increase the overall sensitivity by orders of magnitude**.

Here we introduce a new instrument, the **PTR-TOF 1000 ultra**, which is a compact PTR-TOFMS equipped with an **ion funnel**. Its high sensitivity makes the PTR-TOF 1000 ultra a perfect tool for applications where complex matrices have to be **rapidly analyzed** for example in quality assurance in the food and flavor industry.

One such application is the real-time quantification of **impurities in food-grade carbon dioxide (CO₂)**, which we present here as a practical example.

Impurities in Food-grade CO₂


The modern food and beverage industry is not imaginable without food-grade carbon dioxide. Its many uses include carbonation of beverages, and as packaging and cooling gas. Therefore, we investigated **volatile organic impurities** emerging from a **food-grade certified** (European standard E 290) **carbon dioxide cylinder** with a PTR-TOFMS instrument.

In a first step we compared the high-resolution mass-spectra from the analysis of the CO₂ cylinder a) unfiltered and b) filtered by a charcoal filter (see experimental setup, **left**). The filtered gas serves as a reference for the background of a VOC-free CO₂ matrix. **Figure 4** illustrates that most compounds are suppressed during filtration and therefore constitute **organic impurities, most with mass-to-charge ratios higher than m/z 150 and concentrations in the range of 100 pptv to 10 ppbv**.

In a next step we monitored the concentrations of several main impurities over time (**figure 5**). Immediately after opening the cylinder.


Proton-Transfer-Reaction (PTR)

- ▶ "Soft" ionization of VOCs
 - low fragmentation
 - important in complex matrix
- ▶ Accuracy < 5 %, precision < 2 %
- ▶ Real-time < 100 ms
- ▶ High sensitivity, LOD < 1 ppt



Novel ION BOOSTER


- ▶ Improved ion transfer
 - up to x10 more sensitivity
 - Ion-funnel between PTR ionization and TOF
- ▶ Upgrade of existing PTR-MS
 - Modular design
- ▶ Perfect integration
 - PTR, funnel, TOF, all manufactured in-house



The ion-funnel, consisting of a series of lenses with decreasing apertures and applied RF voltages, focuses the ion trajectories.

Time-of-Flight (TOF)

- ▶ Full spectrum in a split-second
 - Separation of isobars
 - Identification of chemical composition
- ▶ High mass resolution
 - Separation of isobars
 - Identification of chemical composition



▶ Lower LOD (< 5 ppt)




Figure 1: LoDs at m/z 147 for PTR-TOF 1000 and PTR-TOF 1000 ultra for integration times between 0.1 and 60s. LoDs are calculated using the 3σ (standard deviation) method.

▶ Higher sensitivity across entire mass scale




Figure 2: Comparison of the sensitivities for PTR-TOF 1000 and PTR-TOF 1000 ultra. The sensitivities are obtained utilizing compounds from certified gas standards.

▶ Better S/N, faster measurement




Figure 3: Considerably less integration time (s) needed to achieve a relative error of 10% (15%) for PTR-TOF 1000 ultra than for PTR-TOF 1000, while analyzing a concentration of 100 pptv.

▶ Compact dimensions



PTR-TOF 1000 ultra




Figure 4: Comparing mass spectra of filtered (=background) and unfiltered CO₂. Several impurities in food-grade CO₂ are visible.




Figure 5: Concentration vs. time. Impurities show different dynamic, which is attributed to their vapor pressure in the cylinder.

We observed that only one compound (m/z 69.034) reached a constant concentration. Other compounds showed a dynamic behavior, with maxima at different times. We attribute this behavior to the outgassing head-space concentration inside the gas cylinder, which depends on the individual vapor pressure and outgassing dynamics from the liquid CO₂.

Consequently, the time between opening the CO₂ cylinder and introduction of the gas into food can have an influence on the level of impurities measure in the QA as well as in the final product.

These results emphasizes the great advantage of real-time analysis in food & flavor quality assurance.